PLD-010

Personal Locating Device

User Manual

The PLD-010 is a **Personal Locating Device** worn on the wrist and used to locate persons with disabilities (such as Alzheimer's or Autism) that may wander or become lost. The device enables tracking and locating of the patient by trained personnel.

Installation Instructions:

Install Battery with + side into the battery lid. At this point, it is advisable to also inspect the condition of the O-ring on the lid. If the O-ring is cracked, torn, or otherwise damaged, replace O-ring before installing the battery.

Thread lid (with battery installed) into bottom of transmitter body. Tighten with provided tool. PLD-010 will begin transmitting as soon as battery is installed.

With the battery installed. Test that the device is transmitting by listening to it on a receiver, or by using a transmitter tester. **Note:** *Checking that the device is transmitting should be done on a daily basis by a trained care giver.*

Install medical bracelet through the slots on the PLD's body.

Install device onto wrist. Either left or right arm may be used. Battery should be replaced every 30 days.

Battery:

The PLD-010 is powered by a CR2450 Lithium battery. Battery life is 30+ days with a new battery installed.

Transmitting Range

The PLD-010 is a LPRS transmitting device. Typical range is 2 miles in a ground search, or 10+ miles when searching while airborne (from Helicopter or Airplane). Actual range may vary depending on many factors such as height of transmitter, height of receiver, type of terrain and obstructions, and radio noise or interference. Higher ground always gives better range.

Regulatory Statements

For Users in the United States

This Device may not interfere with TV reception or Federal Government radar

This device complies with part 15 of the FCC Rules. Operation is subject to the following 2 conditions

(1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

For users in Canada

Operation is subject to the following two conditions: (1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.

The user/purchaser of this device should take note that as digital TV broadcasting stations are introduced in the geographical area, the medical telemetry links may have to be shifted to operate in other (unused) TV channels. The user/purchaser should ensure that the radio device can adapt to this mode of operation.

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If TV channel 13 is used in the area, the installer shall reduce or adjust the radio frequency radiated power so that nearby TV channel 13 receivers do not receive radio interference from the system installed. Suggestions: A test with a TV receiver equipped with a "rabbit-ear antenna" and tuned to channel 13 should be conducted at the perimeter of the user's intended coverage area and should not overlap other users' areas without the latter's consent. If this does not solve the problem, a channel near the 217 MHz edge and not near 216 MHz should be tried.

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